

## Biotix<sup>®</sup> uTIP<sup>™</sup> Ergonomic Pipette Tips

Compatible with universal pipettes including Gilson<sup>®</sup>, Eppendorf<sup>®</sup>, Biohit<sup>®</sup>, and more!

### Biotix<sup>®</sup> uTIP<sup>™</sup> Pipette Tips, Filtered, Low-retention, 960 tips per pack

Volume	2 µL	10 µL	10 µL XL	20 µL	100 µL	200 µL	300 µL	1000 µL	1250 µL
Sterilized	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Racked	732-2384	732-1429	732-1430	732-1432	732-1433	732-3397	732-1434	732-3398	732-2383

### Biotix<sup>®</sup> uTIP<sup>™</sup> Pipette Tips, Non-Filter, Low-retention, 960 tips per pack

Volume	10 µL		10 µL XL		200 µL		250 µL		300 µL		1000 µL		1250 µL	
Sterilized	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes
Racked	613-2813	613-1376	613-2814	613-1377	613-2815	613-1378	613-5858	613-5857	613-2816	613-1379	613-5863	613-5862	613-2817	613-2436
CleanPak	613-5616	613-5621	613-5615	613-5620	613-5614	613-5619	613-5860	613-5859	613-5613	613-5618	613-5865	613-5864	613-5612	613-5617
Bulk	613-5622		613-5623		613-5624		613-5861		613-5625		613-5866		613-5647	


#### Secure Seal

##### FlexFit<sup>®</sup>

##### Flexible Proximal End

- Gives a consistently secure seal with the widest range of industry leading pipettes
- Requires less insertion force than competitors' to make a secure seal, reducing the risk of repetitive stress injury
- Tips don't fall off your multichannel pipette
- Makes every pipette more ergonomic

## Tip Innovation by Biotix<sup>®</sup>



#### Clean Delivery

##### X-Resin<sup>®</sup>

##### Low Retention Surface

- Maximizes sample recovery
- Improves accuracy across your experiments

##### Blade<sup>®</sup>

##### Minimizes Surface Tension

- Eliminates hanging droplet formation
- Saves time by reducing the need for tip touch off
- Increases accuracy and precision